

Instrument	Parameter	Low	Normal	High
EDAN	WBC × 10 <sup>9</sup> /L	<b>3.43</b> ± 0.60	<b>8.38</b> ± 1.00	<b>18.73</b> ± 2.50
	NEU# × 10 <sup>9</sup> /L	1.65 ± 1.0	4.48 ± 3.80	11.59 ± 7.60
	NEU%	47.8 ± 10.0	53.4 ± 9.0	61.9 ± 8.0
H60	LYM# × 10 <sup>9</sup> /L	1.2 ± 1.0	2.3 ± 1.5	3.4 ± 2.0
H66	LYM	35.2 ± 10.0	27.2 ± 9.0	17.9 ± 8.0
H68	MON# × 10 <sup>9</sup> /L	0.33 ± 0.20	0.76 ± 0.50	1.52 ± 1.10
H69	MON	9.8 ± 4.0	9.1 ± 5.0	8.1 ± 6.0
H60S	EOS# × 10 <sup>9</sup> /L	0.24 ± 0.20	0.83 ± 0.50	2.24 ± 1.30
H66S	EOS	7.2 ± 5.0	9.9 ± 6.0	11.9 ± 7.0
H68S	BAS# × 10 <sup>9</sup> /L	1.0 ± 1.0	1.0 ± 1.0	1.0 ± 1.0
H69S	BAS	1.0 ± 1.0	1.0 ± 1.0	1.0 ± 1.0
	RBC × 10 <sup>12</sup> /L	<b>2.34</b> ± 0.30	<b>4.60</b> ± 0.40	<b>5.25</b> ± 0.50
	HGB g/L	<b>64</b> ± 6	<b>135</b> ± 8	<b>169</b> ± 8
	HCT %	18.6 ± 3.0	40.9 ± 4.0	50.5 ± 5.0
	MCV fL	<b>79.7</b> ± 5.0	<b>88.9</b> ± 5.0	<b>96.2</b> ± 6.0
	MCH pg	27.1 ± 2.5	29.4 ± 2.5	32.1 ± 2.5
	MCHC g/L	341 ± 30	330 ± 30	334 ± 30
	RDW-CV %	15.9 ± 3.0	14.8 ± 3.0	14.1 ± 3.0
	RDW-SD fL	47 ± 10	48 ± 10	50 ± 12
	PLT × 10 <sup>9</sup> /L	<b>53</b> ± 30	<b>270</b> ± 45	<b>525</b> ± 70
	MPV fL	10.7 ± 3.0	10.0 ± 3.0	9.5 ± 3.0
	PDW fL	13.3 ± 4.0	12.1 ± 4.0	11.5 ± 4.0
	PCT %	0.06 ± 0.06	0.27 ± 0.20	0.50 ± 0.20
	PLCR %	30.7 ± 10.0	25.4 ± 10.0	21.4 ± 10.0
	PLCC × 10 <sup>9</sup> /L	16.4 ± 15	68.5 ± 20	112.1 ± 30

## 【NOTE】

- 1) The controls should be stored in the refrigerator (2~8°C). After opening, it will keep stable for 14 days when it is stored airtight at 2~8°C.
- 2) Please equilibrate the controls to room temperature (15~30°C) before using it.
- 3) Controls must be well mixed before use. Please mix gently to avoid cells rupturing or generating bubbles.
- 4) After using, put the controls back into the refrigerator (2~8°C) to prevent contamination and evaporation.